


Product Datasheet

TASP1 antibody (orb1323806)

Description	Mouse monoclonal antibody to TASP1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human TASP1 (NP_060184) produced in HEK293T cell.
Target	TASP1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	44.3 kDa
Uniprot ID	Q9H6P5
NCBI	NM_017714
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...