


Product Datasheet

MSI1 antibody (orb1323808)

Description	Mouse monoclonal antibody to MSI1
Species/Host	Mouse
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human MSI1 (NP_002433) produced in HEK293T cell.
Target	MSI1
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	38.9 kDa
Uniprot ID	O43347
NCBI	NM_002442
Dilution Range	WB 1:1000~2000, IHC 1:150, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...