


Product Datasheet

SNTG1 antibody (orb1323837)

Description	Mouse monoclonal antibody to SNTG1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Full length human recombinant protein of human SNTG1 (NP_061840) produced in HEK293T cell.
Target	SNTG1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	57.8 kDa
Uniprot ID	Q9NSN8
NCBI	NM_018967
Dilution Range	WB 1:2000, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...