


Product Datasheet

Gephyrin (GPHN) antibody (orb1323843)

Description	Mouse monoclonal antibody to Gephyrin (GPHN)
Species/Host	Mouse
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human GPHN (NP_065857) produced in HEK293T cell.
Target	GPHN
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	83.3 kDa
Uniprot ID	Q9NQX3
NCBI	NM_020806
Dilution Range	WB 1:500~2000, IHC 1:150, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...