


Product Datasheet

Josephin 1 (JOSD1) antibody (orb1323850)

Description	Mouse monoclonal antibody to Josephin 1 (JOSD1)
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human JOSD1 (NP_055691) produced in HEK293T cell.
Target	JOSD1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	23 kDa
Uniprot ID	Q15040
NCBI	NM_014876
Dilution Range	WB 1:500, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...