


Product Datasheet

ERK5 (MAPK7) antibody (orb1323905)

Description	Mouse monoclonal antibody to ERK5 (MAPK7)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Full length human recombinant protein of human MAPK7 (NP_002740) produced in HEK293T cell.
Target	MAPK7
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	88.2 kDa
Uniprot ID	Q13164
NCBI	NM_002749
Dilution Range	WB 1:2000, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...