


Product Datasheet

GRHPR antibody (orb1323961)

Description	Mouse monoclonal antibody to GRHPR
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Full length human recombinant protein of human GRHPR (NP_036335) produced in HEK293T cell.
Target	GRHPR
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	35.5 kDa
Uniprot ID	Q9UBQ7
NCBI	NM_012203
Dilution Range	WB 1:2000, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...