


Product Datasheet

PYCR3 antibody (orb1324011)

Description	Mouse monoclonal antibody to PYCR3
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human PYCRL(NP_075566) produced in HEK293T cell.
Target	PYCR3
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	28.5 kDa
Uniprot ID	Q53H96
NCBI	NM_023078
Dilution Range	WB: 1:200 - 1:1000~2000, IHC 1:150, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...