


Product Datasheet

RABL2A antibody (orb1324021)

Description	Mouse monoclonal antibody to RABL2A
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human RABL2A(NP_009013) produced in HEK293T cell.
Target	RABL2A
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	25.9 kDa
Uniprot ID	Q9UBK7
NCBI	NM_007082
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...