


Product Datasheet

HDAC1 antibody (orb1324039)

Description	Mouse monoclonal antibody to HDAC1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Full length human recombinant protein of human HDAC1 (NP_004955) produced in HEK293T cell.
Target	HDAC1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	54.9 kDa
Uniprot ID	Q13547
NCBI	NM_004964
Dilution Range	WB 1:2000, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...