


Product Datasheet

PYCR2 antibody (orb1324040)

Description	Mouse monoclonal antibody to PYCR2
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	Full length human recombinant protein of human PYCR2 (NP_037460) produced in HEK293T cell.
Target	PYCR2
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	33.5 kDa
Uniprot ID	Q96C36
NCBI	NM_013328
Dilution Range	WB 1:2000, IHC 1:150, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...