


Product Datasheet

LRATD2 antibody (orb1324042)

Description	Mouse monoclonal antibody to LRATD2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Full length human recombinant protein of human FAM84B (NP_777571) produced in HEK293T cell.
Target	LRATD2
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	34.3 kDa
Uniprot ID	Q96KN1
NCBI	NM_174911
Dilution Range	WB 1:2000, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...