


Product Datasheet

PECR antibody (orb1324048)

Description	Mouse monoclonal antibody to PECR
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Full length human recombinant protein of human PECR (NP_060911) produced in HEK293T cell.
Target	PECR
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	32.4 kDa
Uniprot ID	Q9BY49
NCBI	NM_018441
Dilution Range	WB 1:500~2000, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...