


Product Datasheet

SPR antibody (orb1324050)

Description	Mouse monoclonal antibody to SPR
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human SPR (NP_003115) produced in HEK293T cell.
Target	SPR
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	27.9 kDa
Uniprot ID	P35270
NCBI	NM_003124
Dilution Range	WB 1:500~2000, IHC 1:150, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...