

Product Datasheet

Bestrophin 3 (BEST3) antibody (orb1324053)

Description Mouse monoclonal antibody to Bestrophin 3 (BEST3)

Species/Host Mouse

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications FC, IF, WB

Immunogen Full length human recombinant protein of human BEST3 (NP_689652) produced in HEK293T cell.

Target BEST3

Note For research use only

Isotype IgG1

Clonality Monoclonal


MW 50.8 kDa

Uniprot ID [Q8N1M1](#)

NCBI [NM_152439](#)

Dilution Range WB 1:2000, IF 1:100, FLOW 1:100

Expiration Date 12 months from date of receipt.



immunofluorescent
staining of COS7
cells...