


Product Datasheet

ZIM2 antibody (orb1324057)

Description	Mouse monoclonal antibody to ZIM2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-150 and 428-527 of human ZIM2 (NP_056178) produced in E.coli.
Target	ZIM2
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	61 kDa
Uniprot ID	Q9NZV7
NCBI	NM_015363
Dilution Range	WB 1:2000, IHC 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...