


Product Datasheet

CALCOCO2 antibody (orb1324060)

Description	Mouse monoclonal antibody to CALCOCO2
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Full length human recombinant protein of human CALCOCO2 (NP_005822) produced in HEK293T cell.
Target	CALCOCO2
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	52.1 kDa
Uniprot ID	Q13137
NCBI	NM_005831
Dilution Range	WB 1:1000~2000, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...