


Product Datasheet

PPWD1 antibody (orb1324083)

Description	Mouse monoclonal antibody to PPWD1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human PPWD1 (NP_056157) produced in HEK293T cell.
Target	PPWD1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	73.4 kDa
Uniprot ID	Q96BP3
NCBI	NM_015342
Dilution Range	WB 1:200~500, IF 1:100, FLOW 1:100, IHC: 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...