


Product Datasheet

BIN3 antibody (orb1324111)

Description	Mouse monoclonal antibody to BIN3
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Full length human recombinant protein of human BIN3 (NP_061158) produced in HEK293T cell.
Target	BIN3
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	29.5 kDa
Uniprot ID	Q9NQY0
NCBI	NM_018688
Dilution Range	WB 1:500, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...