


Product Datasheet

Ribonuclease Inhibitor (RNH1) antibody (orb1324142)

Description	Mouse monoclonal antibody to Ribonuclease Inhibitor (RNH1)
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	Full length human recombinant protein of human RNH1 (NP_002930) produced in HEK293T cell.
Target	RNH1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	49.8 kDa
Uniprot ID	P13489
NCBI	NM_002939
Dilution Range	WB 1:500~2000, FLOW 1:100, IHC: 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...