


Product Datasheet

AIMP2 antibody (orb13242)

Description	AIMP2 antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ICC
Immunogen	KLH conjugated synthetic peptide derived from human AIMP2 (201-300/320 aa)
Target	AIMP2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	36 kDa
Uniprot ID	Q13155
Dilution Range	ICC=1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of rat brai...