


Product Datasheet

TCP11L2 antibody (orb1324207)

Description	Mouse monoclonal antibody to TCP11L2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Full length human recombinant protein of human TCP11L2 (NP_689985) produced in HEK293T cell.
Target	TCP11L2
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	57.9 kDa
Uniprot ID	Q8N4U5
NCBI	NM_152772
Dilution Range	WB 1:500~2000, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...