


Product Datasheet

ERCC1 antibody (orb1324237)

Description	Mouse monoclonal antibody to ERCC1
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Full length human recombinant protein of human ERCC1 (NP_973730) produced in HEK293T cell.
Target	ERCC1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	32.4 kDa
Uniprot ID	P07992
NCBI	NM_001983
Dilution Range	WB 1:2000, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...