


Product Datasheet

IVD antibody (orb1324261)

Description	Mouse monoclonal antibody to IVD
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human IVD (NP_002612) produced in HEK293T cell.
Target	IVD
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	46.5 kDa
Uniprot ID	P26440
NCBI	NM_002225
Dilution Range	WB 1:500, IHC 1:150, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...