


Product Datasheet

DNAJA2 antibody (orb1324293)

Description	Mouse monoclonal antibody to DNAJA2
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Full length human recombinant protein of human DNAJA2 (NP_005871) produced in HEK293T cell.
Target	DNAJA2
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	45.6 kDa
Uniprot ID	O60884
NCBI	NM_005880
Dilution Range	WB 1:500~2000, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...