


Product Datasheet

CBWD1 antibody (orb1324366)

Description	Mouse monoclonal antibody to CBWD1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Full length human recombinant protein of human CBWD1 (NP_060961) produced in HEK293T cell.
Target	CBWD1
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	43.9 kDa
Uniprot ID	Q9BRT8
NCBI	NM_018491
Dilution Range	WB 1:2000, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...