


Product Datasheet

KLHL2 antibody (orb1324369)

Description	Mouse monoclonal antibody to KLHL2
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-100 and 494-593 of human KLHL2 (NP_009177) produced in HEK293T cell.
Target	KLHL2
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	65.8 kDa
Uniprot ID	O95198
NCBI	NM_007246
Dilution Range	WB 1:500, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...