


Product Datasheet

RNF86 (TRIM2) antibody (orb1324421)

Description	Mouse monoclonal antibody to RNF86 (TRIM2)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-100 and 645-744 of human TRIM2 (NP_056086) produced in HEK293T cell.
Target	TRIM2
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	81.3 kDa
Uniprot ID	Q9C040
NCBI	NM_015271
Dilution Range	WB 1:2000, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...