


Product Datasheet

UCK (UCK1) antibody (orb1324438)

Description	Mouse monoclonal antibody to UCK (UCK1)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human UCK1 (NP_113620) produced in HEK293T cell.
Target	UCK1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	31.3 kDa
Uniprot ID	Q9HA47
NCBI	NM_031432
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...