


Product Datasheet

PAICS antibody (orb1324502)

Description	Mouse monoclonal antibody to PAICS
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Full length human recombinant protein of human PAICS (NP_006443) produced in HEK293T cell.
Target	PAICS
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	46.9 kDa
Uniprot ID	P22234
NCBI	NM_006452
Dilution Range	WB 1:2000, IF 1:100, Flow: 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...