


Product Datasheet

mu Crystallin (CRYM) antibody (orb1324535)

Description	Mouse monoclonal antibody to mu Crystallin (CRYM)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human CRYM (NP_001879) produced in HEK293T cell.
Target	CRYM
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	33.6 kDa
Uniprot ID	Q14894
NCBI	NM_001888
Dilution Range	WB: 1:200 - 1:1000, IHC 1:300, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...