


Product Datasheet

Thimet Oligopeptidase (THOP1) antibody (orb1324558)

Description	Mouse monoclonal antibody to Thimet Oligopeptidase (THOP1)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human THOP1 (NP_0003240) produced in HEK293T cell.
Target	THOP1
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	78.7 kDa
Uniprot ID	P52888
NCBI	NM_003249
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...