


Product Datasheet

PRTG antibody (orb1324573)

Description	Mouse monoclonal antibody to PRTG
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Full length human recombinant protein of human PRTG (NP_776175) produced in HEK293T cell.
Target	PRTG
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	126.9 kDa
Uniprot ID	Q2VWP7
NCBI	NM_173814
Dilution Range	WB 1:2000, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...