


Product Datasheet

PEX5 antibody (orb1324577)

Description	Mouse monoclonal antibody to PEX5
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Full length human recombinant protein of human PEX5 (NP_000310) produced in HEK293T cell.
Target	PEX5
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	69.7 kDa
Uniprot ID	P50542
NCBI	NM_000319
Dilution Range	WB 1:2000, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...