


Product Datasheet

COX6B2 antibody (orb1324589)

Description	Mouse monoclonal antibody to COX6B2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Full length human recombinant protein of human COX6B2 (NP_653214) produced in HEK293T cell.
Target	COX6B2
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	10.3 kDa
Uniprot ID	Q6YFQ2
NCBI	NM_144613
Dilution Range	WB 1:500, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...