


Product Datasheet

TTC32 antibody (orb1324599)

Description	Mouse monoclonal antibody to TTC32
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human TTC32 (NP_001008238) produced in HEK293T cell.
Target	TTC32
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	17.1 kDa
Uniprot ID	Q5I0X7
NCBI	NM_001008237
Dilution Range	WB 1:2000, IHC 1:50, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...