


Product Datasheet

HID1 antibody (orb1324651)

Description	Mouse monoclonal antibody to HID1
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human C17orf28 (NP_085133) produced in HEK293T cell.
Target	HID1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	88.6 kDa
Uniprot ID	Q8IV36
NCBI	NM_030630
Dilution Range	WB 1:1000~2000, IHC 1:50, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...