


Product Datasheet

C9orf41 (CARNMT1) antibody (orb1324654)

Description	Mouse monoclonal antibody to C9orf41 (CARNMT1)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human C9orf41 (NP_689663) produced in HEK293T cell.
Target	CARNMT1
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	47.0 kDa
Uniprot ID	Q8N4J0
NCBI	NM_152420
Dilution Range	WB 1:1000~2000, IHC 1:50, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...