


Product Datasheet

CYB5R3 antibody (orb1324666)

Description	Mouse monoclonal antibody to CYB5R3
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human CYB5R3(NP_000389) produced in HEK293T cell.
Target	CYB5R3
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	34.1 kDa
Uniprot ID	P00387
NCBI	NM_000398
Dilution Range	WB 1:500~2000, IHC 1:150, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...