


Product Datasheet

TACC3 antibody (orb1324700)

Description	Mouse monoclonal antibody to TACC3
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human TACC3(NP_006333) produced in HEK293T cell.
Target	TACC3
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	90.2 kDa
Uniprot ID	Q9Y6A5
NCBI	NM_006342
Dilution Range	WB 1:2000, IHC 1:50, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...