

Product Datasheet

Cytochrome P450 1A2 (CYP1A2) antibody (orb1324794)

Description Mouse monoclonal antibody to Cytochrome P450 1A2 (CYP1A2)

Species/Host Mouse

Reactivity Human

Conjugation Unconjugated

Tested Applications IF, WB

Immunogen Full length human recombinant protein of human CYP1A2 (NP_000752) produced in HEK293T cell.

Target CYP1A2

Note For research use only

Isotype IgG1

Clonality Monoclonal


MW 58.2 kDa

Uniprot ID [P05177](#)

NCBI [NM_000761](#)

Dilution Range WB 1:500, IF 1:100

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...