


Product Datasheet

TUBA8 antibody (orb1324826)

Description	Mouse monoclonal antibody to TUBA8
Species/Host	Mouse
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Full length human recombinant protein of human TUBA8 (NP_061816) produced in HEK293T cell.
Target	TUBA8
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	49.9 kDa
Uniprot ID	Q9NY65
NCBI	NM_018943
Dilution Range	WB 1:200~1000, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...