


Product Datasheet

ALDH3A1 antibody (orb1324841)

Description	Mouse monoclonal antibody to ALDH3A1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human ALDH3A1(NP_000682) produced in HEK293T cell.
Target	ALDH3A1
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	50.2 kDa
Uniprot ID	P30838
NCBI	NM_000691
Dilution Range	WB 1:2000, IHC 1:50, IF 1:100, Flow 1:100
Expiration Date	12 months from date of receipt.



Flow
cytometric
HEK293T
cells
transfected...