


Product Datasheet

PIM2 antibody (orb1324882)

Description	Mouse monoclonal antibody to PIM2
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	Full length human recombinant protein of human PIM2 (NP_006866) produced in HEK293T cell.
Target	PIM2
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	34.2 kDa
Uniprot ID	Q9P1W9
NCBI	NM_006875
Dilution Range	WB 1:2000, IHC 1:50, Flow 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...