


Product Datasheet

PTP epsilon (PTPRE) antibody (orb1324895)

Description	Mouse monoclonal antibody to PTP epsilon (PTPRE)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human PTPRE (NP_006495) produced in HEK293T cell.
Target	PTPRE
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	80.5 kDa
Uniprot ID	P23469
NCBI	NM_006504
Dilution Range	WB 1:2000, IHC 1:50, IF 1:100, FLOW 1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...