


Product Datasheet

SIL1 antibody (orb1324903)

Description	Mouse monoclonal antibody to SIL1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human SIL1 (NP_071909) produced in HEK293T cell.
Target	SIL1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	52.1 kDa
Uniprot ID	Q9H173
NCBI	NM_022464
Dilution Range	WB 1:2000, IHC 1:50, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...