


Product Datasheet

Amyloid Precursor Protein (APP) antibody (orb1324939)

Description	Mouse monoclonal antibody to Amyloid Precursor Protein (APP)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Full length human recombinant protein of human APP (NP_000475) produced in HEK293T cell.
Target	APP
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	87.0 kDa
Uniprot ID	P05067
NCBI	NM_000484
Dilution Range	WB 1:2000, IF 1:100, Flow 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...