


Product Datasheet

MCL1 antibody (orb1324961)

Description	Mouse monoclonal antibody to MCL1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human MCL1 (NP_068779) produced in HEK293T cell.
Target	MCL1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	37.3 kDa
Uniprot ID	Q07820
NCBI	NM_021960
Dilution Range	WB 1:2000, IHC 1:50, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...