


Product Datasheet

RALBP1 antibody (orb1324991)

Description	Mouse monoclonal antibody to RALBP1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, IP, WB
Immunogen	Full length human recombinant protein of human RALBP1 (NP_006779) produced in HEK293T cell.
Target	RALBP1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	76.0 kDa
Uniprot ID	Q15311
NCBI	NM_006788
Dilution Range	WB 1:2000, IHC 1:50, IF 1:100, IP: 4ug/mL
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...