


Product Datasheet

XRCC1 antibody (orb1325048)

Description	Mouse monoclonal antibody to XRCC1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human XRCC1 (NP_006288) produced in HEK293T cell.
Target	XRCC1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	69.3 kDa
Uniprot ID	P18887
NCBI	NM_006297
Dilution Range	WB 1:500~2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...